

Branson™ Sonifier Accessories

SFX 150 Microtips

Overview

Microtips for the SFX Sonifier are available in four sizes to handle volumes up to 150 ml. The microtips are also designed to process samples in the smallest diameter vessels.

The 1/8" microtip is already provided with the SFX Sonifiers.



Part #	Processing Volume	Tip Diameter	Approx. Length (working)
109-122-1066	1 ml & under	3/32" (2.4mm)	4"
109-122-1065	3 -10 ml	1/8" (3.2mm)	4"
109-122-1182	5 -25 ml	3/16" (4.8mm)	4"
109-122-1064	10 - 150 ml	1/4" (6.4mm)	4"
109-116-1566	3-10 ml (4) tips	1/8" (3.2mm)	4"

NOTE:
The 109-116-1566 has (4) 1/8" microtips on a 1" center horn.
These tips are for use with the Branson 150 watt SFX and SLP Sonifier only.

BRANSON™ For more information:
www.Emerson.com/Branson

The Emerson logo is a trademark and service mark of Emerson Electric Co. Brand logotype are registered trademarks of one of the Emerson family of companies. All other marks are the property of their respective owners. © 2024 Emerson Electric Co. All rights reserved.

